


BIOSENSE CONTROLLED EC-OWLS MEASUREMENTS

- **Chronoamperometry** – Controlled potential OWLS measurement
- **Cyclic Voltammetry** – either operated alone or simultaneously with OWLS measurement



EC- OWLS measuring setup

CYCLIC VOLTAMMETRY – as a stand alone measurement controlled by BioSense software

Cyclic Voltammetry Settings

Initial potential [Volts]: OK

Initial potential time [s]: Cancel

Start potential [Volts]:


End potential [Volts]:

Sweep rate [Volts/s]:


Volts per point:

Number of Cycles:

Current range:




- 1 Initial potential
- 2 Initial potential time
- 3 Starting potential at time 2 may be equal to 1 or different
- 4 End potential
- 5 Sweep rate
- 6 Volts per point (set the increment of potential where measurements are made)




SIMULTANEOUS OWLS AND CONTROLLED POTENTIAL EC MEASUREMENT

This technique allows applying an external electric field into the flow cell. The electric field dependence of the adsorption/desorption of proteins or other biomolecules can be studied qualitatively and quantitatively by controlled potential OWLS measurement. The value of the potential can be modified during the measurement and the resulted electric current is measured



SIMULTANEOUS OWLS AND CYCLIC VOLTAMMETRY MEASUREMENT

It is possible to run OWLS measurement and Cyclic Voltammetry measurement simultaneously. The BioSense software displays the measured data of both measurements. The CV measurement parameters can be set the same way as in stand alone CV measurement.



References:

1. Stankowski S., Ramsden J.J., Voltage dependent coupling of light into ITO-covered waveguides, *J. Phys. D: Appl. Phys.* 35 (2002) 299-302
2. Brusatori M.A., Van Tassel, P. R., Biosensing under an applied voltage using optical waveguide lightmode spectroscopy, *Biosensors and Bioelectronics* 18 (2003) 1269-1277
3. Bearinger J. P., Vörös J., Hubbell J. A., Textor M., Electrochemical Optical Waveguide Lightmode Spectroscopy (EC-OWLS): A Pilot Study Using Evanescent-Field Optical Sensing under Voltage Control to Monitor Polycationic Polymer Adsorption Onto Indium Tin Oxide (ITO)-Coated Waveguide Chips, *Biotechnology and Bioengineering*, Vol. 82, No. 4, (2003) 465-473