


ITO QCM SENSOR CHIP

for the SRS QCM 200 system (www.thinksrs.com)



FRONT side (sensing side)

BACK side (electrical contact side)

Product description.

Nominal frequency:	5050 ± 10 kHz
Crystal Cutting:	AT
Material:	Pure Z quartz
Angle:	35° 15' ± 30'
XX' Angle:	0° 00' ±30'
Diameter:	25,4 +0/ -0,05 mm
Electrode on the sensor side:	ITO
Electrode on the contact side:	ITO
Apperance:	Transparent in the visible